


## Dual-View - Head to Head



### LBX-02HTBL-1

26" x 42" Table  
22" T-Legs

- 26" x 42" table size with a 28 ¾" work surface height.
- Wilsonart 4168-60-390 Gray Pampus color - chemical safe work surface.
- Charcoal color rubber T-edge around outside of table edge.
- 1 ¼" Laminated wood core with smooth masonry backer board on bottom.
- Oxford color, tubular steel leg supports.
- ¼-20 thread inserts in bottom of table top for mounting legs.

---

#### Olympus Observation Tubes and Leeds Custom Items Required

Qty	Part Number	Description
1	U-DO105	U-DO3-5; Dual Observation Unit, Front-Back Users W/LEDPointer
1	Crating	Leeds Shipping Crate for Table

---

#### Leeds Custom Products

Phone: 800-444-5333 • 763-546-8575


Fax: 763-546-4369

custom@leedsmicro.com

www.leedsmicro.com

by Leeds

## Dual View - Side to Side



### LBX-02HTBL-1

26" x 42" Table

22" T-Legs

- 26" x 42" table size with a 28 ¾" work surface height.
- Wilsonart 4168-60-390 Gray Pampus color - chemical safe work surface.
- Charcoal color rubber T-edge around outside of table edge.
- 1 ¼" Laminated wood core with smooth masonry backer board on bottom.
- Oxford color, tubular steel leg supports.
- ¼-20 thread inserts in bottom of table top for mounting legs.

---

#### Olympus Observation Tubes and Leeds Custom Items Required

Qty	Part Number	Description
1	U-DO123	U-SD05-5; Side Observation Unit, Side By Side User, W/LED Pointer
1	Crating	Leeds Shipping Crate for Table

---

#### Leeds Custom Products


Phone: 800-444-5333 • 763-546-8575

Fax: 763-546-4369

custom@leedsmicro.com

www.leedsmicro.com

# Four Position Observation Microscope Setup



## LBX-02HTBL-1

26" x 42" Table  
22" T-Legs

- 26" x 42" table size with a 28 ¾" work surface height.
- Wilsonart 4168-60-390 Gray Pampus color - chemical safe work surface.
- Charcoal color rubber T-edge around outside of table edge.
- 1 ¼" Laminated wood core with smooth masonry backer board on bottom.
- Oxford color, tubular steel leg supports.
- ¼-20 thread inserts in bottom of table top for mounting legs.

---

### Olympus Observation Tubes and Leeds Custom Items Required

Qty	Part Number	Description
1	U-DO105	U-DO3-5; Dual Observation Unit, Front-Back Users W/LED Pointer
1	U-DO205	U-MDOB3-5 Multiview Observation Unit
1	U-DO301	U-MDOSV-2-5 Multiview Optical Relay Tube - 2 Position
1	Crating	Leeds Shipping Crate for Table

---

### Leeds Custom Products

Phone: 800-444-5333 • 763-546-8575

Fax: 763-546-4369

custom@leedsmicro.com

www.leedsmicro.com

# Five Position Observation Microscope Setup

## LBX-05HTBL-H

22" x 72" Table

20" T-Legs

28" Cross Support Member



- 22" x 72" table size with a 28 ¾" work surface height.
- Wilsonart 4168-60-390 Gray Pampus color, chemical safe work surface.
- Charcoal color rubber T-edge around outside of table edge.
- 1 ¼" Laminated wood core with smooth masonry backer board on bottom.
- 1" square tube for additional support through center length of the table.
- Oxford color, tubular steel leg supports.
- ¼-20 thread inserts in bottom of table top for mounting legs.

---

### Olympus Observation Tubes and Leeds Custom Items Required

Qty	Part Number	Description
1	U-DO205	U-MDOB3-5 Multiview Observation Unit
2	U-DO301	U-MDOSV-2-5 Multiview Optical Relay Tube - 2 Position
1	Crating	Leeds Shipping Crate for Table

---

### Leeds Custom Products

Phone: 800-444-5333 • 763-546-8575


Fax: 763-546-4369

custom@leedsmicro.com

www.leedsmicro.com

by Leeds

# Six Position Observation Microscope Setup



## LBX-05HTBL-H

22" x 72" Table

20" T-Legs

28" Cross Support Member

- 22" x 72" table size with a 28 ¾" work surface height.
- Wilsonart 4168-60-390 Gray Pampus color, chemical safe work surface.
- Charcoal color rubber T-edge around outside of table edge.
- 1 ¼" Laminated wood core with smooth masonry backer board on bottom.
- 1" square tube for additional support through center length of the table.
- Oxford color, tubular steel leg supports.
- ¼-20 thread inserts in bottom of table top for mounting legs.

---

### Olympus Observation Tubes and Leeds Custom Items Required

Qty	Part Number	Description
1	U-DO205	U-MDOB3-5 Multiview Observation Unit
2	U-DO301	U-MDOSV-2-5 Multiview Optical Relay Tube - 2 Position
1	U-DO105	U-DO3-5 Dual Observation - Front to Back Users
1	Crating	Leeds Shipping Crate for Table

---

### Leeds Custom Products

Phone: 800-444-5333 • 763-546-8575

Fax: 763-546-4369

custom@leedsmicro.com

www.leedsmicro.com

# Ten Position Observation Microscope Setup - H-Layout

## LBX-10HTBL-H

48" x 72" Table  
42" Twin-Legs  
60" Cross Support  
Member



- 48" x 72" table size with a 28 ¾" work surface height.
- Wilsonart 4168-60-390 Gray Pampus color, chemical safe work surface.
- Charcoal color rubber T-edge around outside of table edge.
- 1 ¼" Laminated wood core with smooth masonry backer board on bottom.
- 1" square tube for additional support through center length of the table.
- Oxford color, tubular steel leg supports.
- ¼-20 thread inserts in bottom of table top for mounting legs.

---

### Olympus Observation Tubes and Leeds Custom Items Required

Qty	Part Number	Description
1	U-DO212	U-MDO10B3-5 Multiview Observation Unit for 10 Position
4	U-DO301	U-MDOSV-2-5 Multiview Optical Relay Tube - 2 Position
1	U-DO315	U-MDO10R3-5 Rear Projection Unit for 10 Position
1	Crating	Leeds Shipping Crate for Table

---

### Leeds Custom Products

Phone: 800-444-5333 • 763-546-8575

Fax: 763-546-4369

custom@leedsmicro.com

www.leedsmicro.com

# Ten Position Observation Microscope Setup - Straight Layout

## LBX-10HTBL-B2

- 30" x 94" Table
- 20" T-Legs
- 4" Single Post Leg
- 60" Cross Support Member



- 30" x 94" table size with a 28 ¾" work surface height.
- Wilsonart 4168-60-390 Gray Pampus color, chemical safe work surface.
- Charcoal color rubber T-edge around outside of table edge.
- 1 ¼" Laminated wood core with smooth masonry backer board on bottom.
- 1" square tube for additional support through center length of the table.
- Oxford color, tubular steel leg supports.
- ¼-20 thread inserts in bottom of table top for mounting legs.

---

### Olympus Observation Tubes and Leeds Custom Items Required

Qty	Part Number	Description
1	U-DO205	U-MDOB3-5 Multiview Observation Unit
2	U-DO301	U-MDOSV-2-5 Multiview Optical Relay Tube - 2 Position
5	U-DO105	U-DO3-5 Dual Observation. - Front to Back Users
2	LBX-B-0151	Telescope Support Base
2	LBX-153	Support Bracket
1	Crating	Leeds Shipping Crate for Table

---

### Leeds Custom Products

Phone: 800-444-5333 • 763-546-8575

Fax: 763-546-4369

custom@leedsmicro.com

www.leedsmicro.com

# 18 Position Observation Microscope Setup - Open Design



## LBX-18HTBL-2PCS

50" x 96" Table  
 22" T-Legs  
 4" Single Post Legs

- 50" x 96" table size with a 28 ¾" work surface height.
- Wilsonart 4168-60-390 Gray Pampus color, chemical safe work surface.
- Charcoal color rubber T-edge around outside of table edge.
- 1 ¼" Laminated wood core with smooth masonry backer board on bottom.
- Steel plate for additional support through center of the table.
- Oxford color, tubular steel leg supports.
- Table ships in two pieces.
- ¼-20 thread inserts in bottom of table top for mounting legs.

---

### Olympus Observation Tubes and Leeds Custom Items Required

Qty	Part Number	Description
8	U-DO301	U-MDOSV-2-5 Multiview Optical Relay Tube - 2 Position
2	U-DO210	U-MDO10B Modified Multi View Body
2	U-DO203	U-SD03-B-SP; Modified U-SD03 FOR 11+ Mult-Head Configuration
1	U-DO315	U-MDO10R3-5; Rear Projection Unit, 10 Head W/ Led Pointer
1	U-DO205	U-MDOB3-5; Multi Viewing Body
1	Crating	Leeds Shipping Crate for Table

---

### Leeds Custom Products

Phone: 800-444-5333 • 763-546-8575

Fax: 763-546-4369

custom@leedsmicro.com

www.leedsmicro.com

# 26" x 38" Motorized Table



LCT-TBLB shown with LCF-CSTRS



The Leeds Trace Evidence Comparison Microscope on the LCT-TBLB.

The **LCT-TBLB** motorized table is designed to adjust to various heights between 26" for a standard desk chair, or 43" for adjustable stools. (The **LCF-CSTRS** adds 3" to the height of the **LCT-TBLB**).

---

#### Order Product Number

LCT-TBLB	26" x 38" Motorized Table
LCF-CSTRS	Flange Mount Caster Set

---