


Leeds' LEDFO cube illuminator is an LED fiber optic light source designed to meet a wide variety of applications.

- Built with no fan or other moving parts, the Leeds' LED light source is vibration free
- Offers constant color temperature at all light intensities
- Estimated 30,000 hours of operation
- Designed with no thermal component in light emitted
- Accommodates most manufacturers fiber optic ring lights and light guides
- Optional post mount allows LEDFO to be attached to a boom stand to help maximize work space
- Low power consumption



Leeds' LEDFO shown with a stereo microscope, Leeds' boom stand, and Leeds' ring light.



LEDFO Kit Options: LEDFO, Leeds ring light, thread adapter, and Ferule. Part #: LRL-410LEDK



Leeds' LEDFO can be mounted on an optional post to maximize work surface space.