



5-Headed Table



10-Headed Table



18-Headed Table

Leeds offers 3 microscope tables with multiple viewing capabilities. The 5-headed, 10-headed, and 18-headed tables are the most commonly purchased.

TABLE SPECIFICATIONS

Work Surface:	Chemical safe work surface
Color:	
Table Surface Color:	Wilsonart gray color,
Table Trim Color:	Charcoal Color rubber T-edge around table perimeter
Table Underside:	1 1/4" Laminated wood core with smooth masonry backer board underneath
Legs:	Oxford color, tubular steel leg supports
Leg Mounts:	1/4-20 thread inserts in bottom of tabletop
Dimension:	
5-Headed Dimensions:	22" x 72" table size with 28 3/4" work surface height
10-Headed Dimensions:	48" x 72" table size with 28 3/4" work surface height, 42" Twin Legs, 60" Cross Support Member
18-Headed Dimensions:	30" x 94" table size with 28 3/4" work surface height, 20" T Legs, 4" Single Post Leg, 60" Cross Support Member.
18-Headed Table Note:	18-headed table requires extra crating due to size of table.

**Contact Leeds sales staff to discuss your table needs:
E-mail sales@leedsmicro.com or visit
www.leedsmicro.com.**