


Firearms & Tool Marks Comparison Microscope Optional Accessories

Coaxial Illumination Kit




Leeds Coaxial Illumination Kit is an ideal solution for specular reflection, especially when looking at granular breach face marks and analyzing firing pin impressions. Leeds Coaxial Illumination kit provides glare-free illumination.




Breach face viewed through
Leeds LCF3 with Coaxial
Illumination



Breach face
viewed through
Leeds LCF3
with Fluorescent
Illumination



Leed Quick Change Kit Add this mounting for easy swapping out of objectives




LCF-CMOS5-HDR CAMERA

Leeds Special HDR Camera optimized for bullets and cartridge cases. It has a unique high dynamic range mode to eliminate glare and more accurately capture what is seen through the eye pieces.



Leeds Ultra-Widefield Objective


.18x objective, largest FOV at 170mm. Continuous 1.35x to 27x parfocal zoom range. WD of 106.5mm. Easily swap out different objectives with the use of Leeds Quick Change kit.



SAMPLE HOLDERS



Magnetized Shell Holder
Sold as a Pair



White Contrast Disk for Universal Holder
Sold as a Pair



Mega Holder



Wax Platten Holder



Ball Pivot Positioner Large




Ball Pivot Vise Holder



Illuminator Reflector
mounted on flexible tube arm,
Sold as a Pair



On-Center Holder



Reflected Document Holder

POST MOUNTS | SHELVING | PIN MOUNTS & EYEPIECE HOLDER

Monitor Post Mount

Post can be attached to the back of the table and can be oriented either to the right, the left or middle of table.



Keyboard & PC Mounts

Easily attach to the microscope table leg (either right or left) and can be adjusted to desired height.



LCF3 Table Extension Kit


Extension table fits the left or right side of the LCF3-station, includes a neoprene pad insert, and Leeds unique lip anti-roll off lip.



LCF 16" SHELF

LCF 16" shelf attached to the LCF3

LEEDS EYEPIECE AND PIN MOUNT HOLDER




Pin Mount Holder Tray, attaches to eyepiece holder (required when adding on).



Eyepiece Holder, attached to the LCF column & supports the pin mount holder.



Set of Pin Mount, Pairs



LCF 12" SHELF


Eyepiece Set



Includes 10x eyepiece, Widefield focusing 10x, & 10x Crosshair.

LIGHTING OPTIONS

LEEDS FLUORESCENT QUAD-TUBE LIGHTS



LCF2 Fluorescent Kit
Cross Bar Mount




LCF2 Fluorescent Light
Stage Mount


LEEDS LED LIGHT STIX



STAGE MOUNT
articulated arm
attachment, stage
mount



COLUMN MOUNT
flexible tube arm,
column mount.
LCF3 ONLY



CROSS BAR MOUNT
flexible tube arm,
column mount

QUAD-TUBE REPLACEMENT




Fluorescent
Quad-tube,
color-matched
pair of
replacement
Bulbs. Sold as
set of 2.



TWIN-TUBE REPLACEMENT




Twin-tube
Fluorescent
replacement
Bulbs, Sold as
set of 2.



Light Guide Positioner


Used to record and duplicate the distance from sample, height, angle of illumination, and rotary position. Useful for reverse lighting technique.



Bifurcated Fiber Optic Light Guide



Diffuser Lens



Focusing Spot Lens



Front mount & fiber optic cable
shown together

=




Front mount attaches to LCF

+



LCF Fiber Optic Cable

MISCELLANEOUS ACCESSORIES




LCF3-SCALE

Upgraded scale to replace previous versions of LCF scales.



LCF-COVER


Dust Cover for LCF



Leeds Calibration Kit

Magnification matching kit

CASTERS | ANTI-VIBRATION FEET



Anti-vibration pads cut 3"x3" to fit under LCF Legs



Leveling Caster, Sold as set of 4



Anti-vibration feet, Sold as set of 4 for LCF or LCF3



High deflection anti-vibration feet, Sold as set of 4 for LCF3