


Fluorescent Quad Lamp


The Leeds Fluorescent Quad Lamp provides bright diffuse illumination. The unique rotating hood allows examiners to control the angle of light and adjust the contrast on a sample. The Quad Lamp can be attached with an articulated arm, allowing end user to adjust and tighten joints as desired. Alternately, the Quad Lamp can be attached to the microscope mounted on obedient gooseneck tubing, with stay in place capabilities. Adjust the tubing with one hand, and the quad lamp will remain in the desired location without needed to be locked in place. These replacement lamps are quality control checked and sorted by their background color, to minimize the large color variance that can occur between individual lamps. Replacement bulbs may be purchased in pairs on our ebay store. Please click on "Ebay Store" link on our website: www.leedsmicro.com.



Fluorescent Quad Lamp mounted on the LCF3 Firearms and Tool Marks Comparison Microscope



Fluorescent Quad Lamp with articulated arm, attached to Focus Mount



Fluorescent Quad Lamp mounted to column with flexible tube arm



Fluorescent Quad Lamp replacement bulbs, sold in a set of 2



**Please contact Leeds sales staff to discuss your lighting needs:
E-mail: sales@leedsmicro.com, or
visit www.leedsmicro.com.**