



SIX-BUTTON KEY PAD CONTROLLER


ALTERNATIVE TO THE LEEDS TOUCHPAD



Leeds Six- Button Key Pad is a tactile controller for the Leeds LCF3 or Discovery Comparison Microscopes to **control the table height, column height, and fluorescent light controls.**

- Enables the examiner to easily and quickly adjust to multiple positions without moving eyes from the microscope eyepiece.
- With only 6 tactile buttons, your brain easily remembers their positions and functions .
- Leeds Six- Button Key Pad measures 5" x 3.75" x 7/8".

To add this to your current LCF3 or Discovery microscope, or opt for it on your next purchase, contact us: sales@leedsmicro.com



The new tactile Six-Button Controller is an alternative to the touchpad on your Leeds LCF3 or Discovery microscope



LEEDS FORENSIC SYSTEMS
17300 MEDINA ROAD, SUITE 600
MINNEAPOLIS, MN 55447
+1.763-546-8575
WWW.LEEDSMICRO.COM