



LEEDS COAXIAL ILLUMINATION KIT



Leeds Coaxial Illumination Kit is an ideal solution for specular (mirror-like) reflection, especially when looking at granular breachface marks and analyzing firing pin impressions. Leeds Coaxial Illumination kit provides glare-free illumination.



Breachface viewed through Leeds LCF3 microscope with Fluorescent Illumination




Breachface viewed through Leeds LCF3 microscope with Coaxial Illumination

Leeds Coaxial illumination kit includes:


- Powerful LED light source with contrast adjustment control knob.
- Leeds unique quick-change feature incorporating a direct electrical contact.
- Can be retro-fitted to any existing Leeds firearm comparison microscope.
- When not in use, entire kit secures to compact storage station, out of the way of microscope main controls.




LEFT: Leeds Coaxial Kit includes all components pictured.



Power Supply/Storage Station mounts to post at back of microscope table, illuminators mount to power supply and 1/4 wave filters attach to dock on outside of illuminators.



To learn more about this and other Leeds products, please contact Leeds sales staff at sales@leedsmicro.com.



 **LEEDS**
FORENSIC SYSTEMS

LEEDS FORENSIC SYSTEMS
17300 MEDINA ROAD, SUITE 600
MINNEAPOLIS, MN 55447
+1.763-546-8575
WWW.LEEDSMICRO.COM