


# Firearms & Tool Marks Comparison Microscope *Optional Accessories*

## Coaxial Illumination Kit --




an ideal solution for specular reflection, (mirror-like) especially when looking at granular breach face marks and analyzing firing pin impressions. This Kit provides glare-free illumination.




Breach face viewed through  
Leeds LCF3 with Coaxial  
Illumination



Breach face  
viewed through  
Leeds LCF3  
with Fluorescent  
Illumination



**Leed Quick Change Kit** Add this mounting  
for easy swapping out of objectives



## LCF-CMOS5-HDR CAMERA

Leeds Special HDR Camera optimized for bullets and cartridge cases. It has a unique high dynamic range mode to eliminate glare and more accurately capture what is seen through the eye pieces.




## Leeds Ultra-Widefield Objective

.18x objective, largest FOV at 170mm. Continuous 1.35x to 27x parfocal zoom range. WD of 106.5mm. Easily swap out different objectives with the use of Leeds Quick Change kit.






## SAMPLE HOLDERS



**Magnetized Shell Holder**  
Sold as a Pair



**White Contrast Disk for Universal Holder**  
Sold as a Pair



**Mega Holder**



**Wax Platten Holder**



**Ball Pivot Positioner Large**




**Ball Pivot Vise Holder**



**Illuminator Reflector**  
mounted on flexible tube arm,  
Sold as a Pair



**On-Center Holder**



**Reflected Document Holder**

## POST MOUNTS | SHELVING | PIN MOUNTS & EYEPIECE HOLDER

### Monitor Post Mount

Post can be attached to the back of the table and can be oriented either to the right, the left or middle of table.



### Keyboard & PC Mounts

Easily attaches to the microscope table leg (either right or left) and can be adjusted to desired height.



### LCF3 Table Extension Kit


Extension table fits the left or right side of the LCF3-station, includes a neoprene pad insert, and Leeds unique lip anti-roll off lip.



**LCF 16" SHELF**

LCF 16" shelf attached to the LCF3

## LEEDS EYEPIECE AND PIN MOUNT HOLDER




**Pin Mount Holder Tray**, attaches to eyepiece holder (required when adding on).



**Eyepiece Holder**, attached to the LCF column & supports the pin mount holder.



**Set of Pin Mount, Pairs**



**LCF 12" SHELF**


### Eyepiece Set



Includes 10x eyepiece, Widefield focusing 10x, & 10x Crosshair.

# LIGHTING OPTIONS

## LEEDS FLUORESCENT QUAD-TUBE LIGHTS



LCF2 Fluorescent Kit  
Cross Bar Mount



LCF2 Fluorescent Light  
Stage Mount

## LEEDS LED LIGHT STIX



**STAGE MOUNT**  
articulated arm  
attachment, stage  
mount



**COLUMN MOUNT**  
flexible tube arm,  
column mount.  
**LCF3 ONLY**



**CROSS BAR MOUNT**  
flexible tube arm,  
column mount

## QUAD-TUBE REPLACEMENT




Fluorescent  
Quad-tube,  
color-matched  
pair of  
replacement  
Bulbs. Sold as  
set of 2.



## TWIN-TUBE REPLACEMENT




Twin-tube  
Fluorescent  
replacement  
Bulbs, Sold as  
set of 2.



## Light Guide Positioner

Used to record and duplicate the distance from sample, height, angle of illumination, and rotary position. Useful for reverse lighting technique.




Bifurcated Fiber Optic Light  
Guide



Diffuser Lens



Focusing Spot Lens




Reflector attaches to station



Front mount attaches to  
LCF column

+




LCF Fiber Optic Cable

=




Front mount & Fiber  
Optic Cable



Bullet sample without Reflector vs. with



## OTHER ACCESSORIES




### LEED3-SCALE

Upgraded scale to replace previous versions of LCF scales.



### LCF-COVER

Dust Cover for LCF



### Leeds Calibration Kit

Magnification matching kit

## CASTERS | ANTI-VIBRATION FEET



Anti-vibration Pads, sold as set of 4, for LCF2 or LCF3



Leveling Caster, Sold as set of 4



Anti-vibration feet, Sold as set of 4 for LCF or LCF3



High deflection anti-vibration feet, Sold as set of 4 for LCF3