


LIGHT GUIDE POSITIONER

ACCESSORY FOR LIGHTING REPEATABILITY





Leeds Light Guide Positioner (LGP) is an optional accessory to add to your comparison microscope configuration for lighting repeatability.

LGP capabilities include:

- Record and duplicate the distance from the sample, the height and angle of illumination, and the rotary position.
- Reverse Lighting -- the examiner can setup lighting from the same length, height, and angle of illumination of the sample on both sides (mirrored illumination of each other.) This feature assures that both samples are illuminated in the exact same manner.
- Ability to rotate the illumination source around the sample 360 degrees.

LGP on both sides of a Leeds comparison microscope:



LEEDS FORENSIC SYSTEMS, INC.
17300 MEDINA ROAD, SUITE 600
MINNEAPOLIS, MN 55447
+1.763-546-8575
WWW.LEEDSMICRO.COM