




ULTRA-WIDEFIELD OBJECTIVE

Leeds Ultra-Widefield Objective offers:

- An ultra-widefield 0.18x objective.
- The largest field of view (FOV) in the industry at 170mm.
- A continuous 1.35x to 27x parfocal zoom range.
- A working distance of 106.5mm.
- Use on Leeds Discovery-Z or LCF3 Comparison Microscopes.



This objective can be used with Leeds optional Quick Change adapter kit -- allowing user to easily and safely switch out large objectives, sliding each one in/out of a top-mounted "docking station."



To learn more about this and other Leeds products, contact Leeds sales staff at sales@leedsmicro.com.



LEEDS FORENSIC SYSTEMS
17300 MEDINA ROAD, SUITE 600
MINNEAPOLIS, MN 55447
+1.763-546-8575
WWW.LEEDSMICRO.COM