



SIX-BUTTON KEY PAD CONTROLLER

ALTERNATIVE TO THE LEEDS TOUCHPAD




Leeds Six-Button Key Pad is a tactile controller for the Leeds LCF3 or Discovery Comparison Microscopes to **control the table height, column height, and fluorescent light controls.**

- Enables the examiner to easily and quickly adjust to multiple positions without moving eyes from the microscope eyepieces.
- With only 6 tactile buttons, your brain easily remembers their positions and functions.
- Leeds Six-Button Key Pad measures 5" x 3.75" x 7/8".



To add this to your current LCF3 or Discovery microscope, or opt for it on your next purchase, contact us: sales@leedsmicro.com



The new tactile Six-Button Controller is an alternative to the touchpad on your Leeds LCF3 or Discovery microscope



LEEDS PRECISION INSTRUMENTS, INS.
DBA LEEDS FORENSIC SYSTEMS
17300 MEDINA ROAD, SUITE 600
MINNEAPOLIS, MN 55447 USA
763-546-8575 // WWW.LEEDSMICRO.COM