




Light Guide Positioner

A Microscope Accessory for Lighting Repeatability


Leeds compact Light Guide Positioner for your comparison microscope will enhance its functionality.

A compact and elegant solution, the LGP provides the ability to:

- **Record and duplicate the distance** from the sample, the height and angle of illumination, and the rotary position.
- **Reverse Lighting** -- the examiner can setup lighting from the same length, height, and angle of illumination and 180° of the sample on both sides (mirrored illumination of each other). This feature assures that both samples are illuminated equally but opposite.
- Ability to rotate the illumination source around the sample 360°.



Leeds Light Guide Positioner on both sides of the comparison microscope



LEEDS PRECISION INSTRUMENTS, INC.
DBA LEEDS FORENSIC SYSTEMS
17300 MEDINA RD, SUITE 600
MINNEAPOLIS MN 55477 USA
763-546-8575 // WWW.LEEDSMICRO.COM