




Ultra-Widefield Objective for Leeds Discovery Comparison Microscopes

5" Door Strike Plate



Leeds Ultra-Widefield Objective, designed for our Discovery and Discovery-Z provides:


- An ultra-widefield 0.18x objective.
- The largest field of view (FOV) in the industry at 203mm.
- A continuous 1.35x to 27x parfocal zoom range.
- A working distance of 106.5mm.



This objective can be used with Leeds optional **Quick Change adapter kit** (shown at left) -- allowing user to easily and safely switch out large objectives, sliding each one in/out of a top-mounted "docking station."



To learn more about this and other Leeds products, contact sales@leedsmicro.com.



LEEDS PRECISION INSTRUMENTS, INC.
DBA LEEDS FORENSIC SYSTEMS
17300 MEDINA ROAD, SUITE 600
MINNEAPOLIS MN 55447 USA
763-546-8575 // WWW.LEEDSMICRO.COM