



# Optional Accessories for Firearms & Tool Marks Comparison Microscopes


## Coaxial Illumination Kit

Leeds Coaxial Illumination Kit is an ideal solution for specular reflection, especially when looking at granular breach face marks and analyzing firing pin impressions. The kit provides glare-free illumination.



## Leeds Quick-Change Mounting Kit.

Add this to quickly & easily change objectives.



## LCF-CMOS5-PKG Camera

Leeds 3"-tall camera is optimized for bullets and cartridge cases. It has a unique **HDR (high dynamic range)** mode to eliminate glare and more accurately capture what is seen through the eyepieces. The **EDF (extended depth of field)** feature enables the entire sample to be in focus at once.




## Leeds Ultra-Widefield Objective

0.18x objective, largest FOV at 170mm. Continuous 1.35x to 27x parfocal zoom range. WD of 106.5mm.






# SAMPLE HOLDERS



**Magnetized Shell Holder**  
Sold as a Pair



**White Contrast Disk for Universal Holder**  
Sold as a Pair



**Mega Holder**



**Wax Platen Holder**




**Ball Pivot Positioner Large**




**Ball Pivot Vise Holder**



**Illuminator Reflector**  
mounted on flexible tube arm,  
Sold as a Pair



**On-Center Holder**



**Reflected Document Holder**

# POST MOUNTS / SHELVING / PIN MOUNTS & EYEPIECE HOLDER



**Articulated Monitor Arm**  
This arm can be positioned just about anywhere for your convenience!

It attaches to a post that is provided with the Leeds table, on either front or back of the table.




**Keyboard & PC Mounts**  
Easily attach to the microscope table leg (either right or left.) Adjustable to desired height.



**LCF3 Table Extension Kit**  
Extension table fits the left or right side of the LCF3-station, includes a neoprene pad insert, and Leeds unique lip anti-roll off lip.

## LEEDS EYEPIECE AND PIN MOUNT HOLDER



Pin Mount Holder Tray, attaches to eyepiece holder (required when adding on).



Keep extra eyepieces handy with our **eyepiece & pin mount holder**, which attaches to the LCF column .



Set of Pin Mounts, Pairs



LCF leg-mounted shelf, available in 12" or 16" .


**Eyepiece Set**  
includes:  
10x eyepiece;  
Widefield  
focusing 10x;  
10x Crosshair.



**Eyepiece Set**

# LIGHTING

## LEEDS FLUORESCENT QUAD-TUBE LIGHTS




LCF2 Fluorescent Kit  
Cross Bar Mount




LCF2 Fluorescent Light  
Stage Mount


## LEEDS LED STIX



**STAGE MOUNT**  
articulated arm  
attachment, stage  
mount




**COLUMN MOUNT**  
flexible tube arm,  
column mount.  
**LCF3 ONLY**



**CROSS BAR MOUNT**  
flexible tube arm,  
column mount

## QUAD-TUBE REPLACEMENT



Fluorescent  
Quad-tube,  
color-matched  
pair of  
replacement  
Bulbs. Sold as  
set of 2.



## TWIN-TUBE REPLACEMENT




Twin-tube  
Fluorescent  
replacement  
Bulbs, Sold as  
set of 2.



## Light Guide Positioner

Used to record and duplicate the distance from sample, height, angle of illumination, and rotary position. Useful for reverse lighting technique.




Bifurcated Fiber Optic Light  
Guide



Diffuser Lens



Focusing Spot Lens



Light up BOTH sides of the bullet specimen with our reflector accessory. Reflector attaches to station, as shown below.



Front mount attaches to  
LCF column




LCF Fiber Optic Cable



Front mount & Fiber  
Optic Cable




## OTHER ACCESSORIES



**LCF3-SCALE** (Upgraded scale to replace previous versions of LCF scales.)




**LCF-COVER**  
Dust Cover for LCF



**Leeds Calibration Kit**  
Magnification matching kit


## CASTERS / ANTI-VIBRATION FEET

Anti-vibration  
Pads, sold as set of  
4, for LCF2 or LCF3



Leveling Caster,  
Sold as set of 4

Anti-vibration feet,  
Sold as set of 4 for  
LCF or LCF3



High deflection  
anti-vibration feet,  
Sold as set of 4 for  
LCF3