








Ultra-Widefield Objective for Leeds LFC3 & LFC2 Comparison Microscopes

Leeds Ultra-Widefield Objective,
designed for our LFC3 & LFC2, offers:

- An ultra-widefield 0.18x objective.
- The largest field of view (FOV) in the industry at **200mm**.
- A continuous 1.08x to 18x parfocal zoom range.
- A **working distance of 106.5mm**.



This objective can be used with Leeds optional **Quick Change Adapter Kit** (shown above) -- allowing user to easily and safely switch out large objectives, sliding each one in/out of a top-mounted "docking station."



LEEDS PRECISION INSTRUMENTS, INC.
DBA LEEDS FORENSIC SYSTEMS
17300 MEDINA ROAD, SUITE 600
MINNEAPOLIS, MN 55447 USA
763-546-8575 // WWW.LEEDSMICRO.COM