








# Ultra-Widefield Objective for Leeds LFC3 & LFC2 Comparison Microscopes

Leeds Ultra-Widefield Objective, designed for our LFC3 & LFC2, offers:

- An ultra-widefield 0.18x objective.
- The largest field of view (FOV) in the industry at 200mm.
- A continuous 1.08x to 18x parfocal zoom range.
- A working distance of 106.5mm.



This objective can be used with Leeds optional **Quick Change adapter kit** (shown above) -- allowing user to easily and safely switch out large objectives, sliding each one in/out of a top-mounted "docking station."



LEEDS PRECISION INSTRUMENTS, INC.  
DBA LEEDS FORENSIC SYSTEMS  
17300 MEDINA ROAD, SUITE 600  
MINNEAPOLIS, MN 55447 USA  
763-546-8575 // [WWW.LEEDSMICRO.COM](http://WWW.LEEDSMICRO.COM)