



Microscope Service & ISO 17025 Accredited Calibration

Leeds provides the only factory-trained service team for Leeds Comparison Microscopes! Our team is uniquely qualified to keep your instruments in the best possible working order.

We also provide ISO17025-compliant calibration certificates to satisfy any requirements of your quality program.

We offer a complete service program that includes: **calibration, field service, and factory repairs.**

- Leeds services *all* makes and models of microscopes.
- Service can be performed on-site or at Leeds Minneapolis facility.
- Nationwide service.
- We carry a full line of replacement bulbs.

Leeds can provide ISO 17205 accredited, NIST-traceable calibrations for:

- Magnification matching
- Camera measurement functions
- Digital LCD stage micrometers
- Linear eyepiece reticles



Leeds is accredited to ISO/IEC 17025:2017 by PJLA Accreditation Services for dimensional calibration of microscopes.



**PERRY JOHNSON LABORATORY
ACCREDITATION, INC.**

Certificate of Accreditation

Perry Johnson Laboratory Accreditation, Inc. has assessed the Laboratory of:

Leeds Precision Instruments, Inc.
17300 Medina Road, Suite 600, Minneapolis, MN 55447

(Hereinafter called the Organization) and hereby declares that Organization is accredited in accordance with the recognized International Standard:

**ISO/IEC 17025:2017
& Meets Requirements of ANSI/NCSL Z540-1-1994**

This accreditation demonstrates technical competence for a defined scope and the operation of a laboratory quality management system (as outlined by the joint ISO-ILAC-IAF Communiqué dated April 2017):

Dimensional Calibration
(As detailed in the supplement)

Accreditation claims for such testing and/or calibration services shall only be made from addresses referenced within this certificate. This Accreditation is granted subject to the system rules governing the Accreditation referred to above, and the Organization hereby covenants with the Accreditation body's duty to observe and comply with the said rules.

For PJLA:

	<small>Initial Accreditation Date:</small> June 29, 2015	<small>Issue Date:</small> October 18, 2021
<small>Tracy Szerszen President</small>	<small>Revision Date:</small> August 22, 2022	<small>Expiration Date:</small> December 31, 2023
	<small>Accreditation No.:</small> 85567	<small>Certificate No.:</small> L21-630-R1

Perry Johnson Laboratory Accreditation, Inc. (PJLA)
755 W. Big Beaver, Suite 1325
Troy, Michigan 48064

The validity of this certificate is maintained through ongoing assessments based on a continuous accreditation cycle. The validity of this certificate should be confirmed through the PJLA website: www.pjlabs.com

Page 1 of 2




**Call or visit our website to easily schedule your microscope service.
763-546-8575 // www.leedsmicro.com // sales@leedsmicro.com**

Like any complex mechanical instrument, a Leeds comparison microscope needs regular service & preventive maintenance.

- Internal components with grease require re-lubrication over time;
- Belts and gears can potentially deteriorate and fail;
- Internal optics can become dusty.

In the absence of regular, comprehensive maintenance, lubrication can dry out on friction surfaces -- which can cause wear & tear on the microscope components. This can lead to mechanical failure. Also, dust accumulation in the optical path will diminish the quality of your image.

Request microscope service on our website:
www.leedsmicro.com/request-service/



Leeds Precision Instruments, Inc.
DBA Leeds Forensic Systems
17300 Medina Road, Suite 600
Minneapolis, MN 55447 USA
763-546-8575 // www.leedsmicro.com