



Microscope Service

Leeds offers service, preventive maintenance, and repair of microscopes.

Why? Over time, lubrication can dry on the friction surfaces of internal microscope components. This can lead to wear & tear, and potential parts failure. Dirt, dust, grease, and immersion oil may settle on the external & internal optical surfaces of your scope, degrading its ability to present a usable image. A regular service program will keep your microscopes in the best possible working condition and extend the life of the instruments.



Leeds comprehensive microscope service includes:

- Cleaning of the internal & external optical path.
- Re-lubrication of components that hold dried out lube -- including eyepiece diopters, stages, focus assemblies, condensers, sample positioners, and stage clips.
- Verification & adjustment of mechanical function.
- Alignment and centration of optical path.
- Cleaning of microscope body.
- Minor on-site repairs that restore functionality but do not require new parts.

Leeds Service Team is an authorized service provider of Olympus microscopes -- and has over 30 years experience servicing and repairing all major microscope brands.



On Leeds imaging systems, the technician may get into the heart of the system to ensure proper working order.

Call 763-546-8575 or visit www.leedsmicro.com to easily schedule your microscope service

Leeds is accredited to ISO/IEC 17025:2017 by PJLA Accreditation Services for dimensional calibration of microscopes.



- Leeds services **all makes and models** of microscopes.
- Service can be performed on-site or at our Minneapolis facility.
- Leeds provides service nationwide.
- We carry a full line of replacement bulbs. See our Ebay store, linked on www.leedsmicro.com.

Ask about our:

- Annual and Multi-year Service Packages
- Equipment Relocation Services
- On-Site Service by trained and certified technicians
- ISO/IEC 17025:2017 Accreditation for calibrations (for applicable products)
- Customer Service Support and troubleshooting
- Software/Hardware updates
- On-Site Training support



Leeds Precision Instruments, Inc.
DBA Leeds Forensic Systems
17300 Medina Road, Suite 600
Minneapolis, MN 55447 USA
763-546-8575 // www.leedsmicro.com