




Universal Holder Contrast Disc Accessory LCF-BDWT



LCF-BDWT Contrast Disc provides a white backdrop to a bullet or cartridge case when loaded in the Leeds Universal Holder and viewed through a Leeds comparison microscope.

The Contrast Disc is secured with a screw and can be swiveled in and out of viewing range as needed.



LCF-BDWT Contrast Disc Features:

- The white background contrast disc attachment for the Leeds Universal Holder creates a consistent white backdrop for any object viewed and captured through a Leeds comparison microscope.
- Sold as a Pair of 2 contrast discs, one disc for each universal holder on a Leeds Comparison Microscope.
- Part #LCF-BDWT. Contact sales@leedsmicro.com