

LCF3 Can Help Get Actionable Results *Faster*



LCF3



The LCF3 comparison microscope was designed incorporating examiner feedback, focusing on image resolution, sample handling, ergonomics, and versatility. These design concepts provide examiners with **features that can increase the rate at which examiners conduct forensic examinations** --- allowing for **faster actionable results**, such as helping move through

case backlog or providing quicker hit confirmations from automated ballistics identification systems like Evofinder®

This powerful workhorse is used by police department firearms labs across the world as an essential investigative tool in firearm investigations.

LCF3 is available via GSA1122.

The LCF3, a reliable and versatile comparison microscope offers:

- Continuous, parfocal and parcentric primary magnification range from 6x-102x, with apochromatically corrected optics, and 14 magnification set points.
- Durable build with easily interchangeable parts.
- Manual system designed around a truly universal holder for fast inspection of evidence.
- Rail system for mounting accessories.
- 5 year warranty.
- All magnification matching and measurements from the microscope are ISO/IEC 17025:2017 calibrated.



LEEDS PRECISION INSTRUMENTS, INC.
17300 MEDINA RD, SUITE 600
MINNEAPOLIS MN 55447 USA
763-546-8575
WWW.LEEDSMICRO.COM