




# Optional Accessories for Firearms & Tool Marks Comparison Microscopes





**Coaxial Illumination Kit** -- an ideal solution for specular reflection, (mirror-like) especially when looking at granular breach face marks and analyzing firing pin impressions. This Kit provides glare-free illumination.



Breach face viewed through Leeds LCF3 with Coaxial Illumination




Breach face viewed through Leeds LCF3 with Fluorescent Illumination



### Leeds Quick Change Kit

Add this mounting kit to easily and quickly swap out objectives



### LCF-CMOS5-HDR CAMERA

Leeds HDR Camera is optimized for the rounded shape of bullets and cartridge cases. It has a unique high dynamic range mode to eliminate glare and more accurately capture what is seen through the eye pieces.


### Leeds Ultra-Widefield Objective

**For Discovery Microscopes:**  
0.18x objective, largest FOV at 170mm. Continuous 1.35x to 27x parfocal zoom range. Working Distance of 106.5mm.


**For LCF-3 Microscopes:**  
0.18x objective, largest FOV at 200mm. Continuous 1.08x to 18x parfocal zoom range. Working Distance of 106.5mm.



# SAMPLE HOLDERS



**Magnetized Shell Holder**  
Sold as a Pair




**White Contrast Disk for Universal Holder** (Sold as a Pair)



**Mega Holder**



**Wax Platten Holder**




**Ball Pivot Positioner with 4" diameter platen**




**Ball Pivot Vise Holder**



**Illuminator Reflector** mounted on flexible tube arm.  
(Sold as a Pair)




**On-Center Holder**



**Reflected Document Holder**

# POST MOUNTS / SHELVING / PIN MOUNTS & EYEPIECE HOLDER



**Monitor Post Mount**

Post can be attached to the back of the table and can be oriented either to the right, the left, or middle of table.



**Keyboard & PC Mounts**

Easily attaches to the microscope table leg (either right or left) and can be adjusted to desired height.



**LCF3 Table Extension Kit**


Extension table fits the left or right side of the LCF3-station, includes a neoprene pad insert, and Leeds unique lip anti-roll off lip.




**LCF 16" SHELF**

LCF 16" shelf attached to the LCF3

## LEEDS EYEPIECE AND PIN MOUNT HOLDER



**Pin Mount Holder Tray**, attaches to eyepiece holder; required when adding on.



**Eyepiece Holder**, attached to the LCF column & supports the pin mount holder.



**Set of Pin Mounts (Pairs)**



**LCF 12" SHELF**


## Eye-piece Set



Includes (a) 10x eyepiece, (b) Widefield focusing 10x, and (c) 10x Crosshair

# LIGHTING OPTIONS

## LEEDS FLUORESCENT QUAD-TUBE LIGHTS




LCF2 Fluorescent Kit  
Cross Bar Mount




LCF2 Fluorescent Light  
Stage Mount

## LEEDS LED LIGHT STIX



**STAGE MOUNT**  
articulated arm  
attachment, stage  
mount




**COLUMN MOUNT**  
flexible tube arm,  
column mount.  
LCF3 ONLY



**CROSS BAR MOUNT**  
flexible tube arm,  
column mount

## QUAD-TUBE REPLACEMENT



Fluorescent  
Quad-tube,  
color-matched  
pair of  
replacement  
Bulbs. Sold as  
set of 2.



## TWIN-TUBE REPLACEMENT




Twin-tube  
Fluorescent  
replacement  
Bulbs, Sold as  
set of 2.



### Light Guide Positioner

Used to record and duplicate the distance from sample, height, angle of illumination, and rotary position. Useful for reverse lighting technique.




Bifurcated Fiber Optic Light  
Guide



Diffuser Lens




Focusing Spot Lens



Reflector attaches to station




Front mount attaches  
to LCF column



LCF Fiber Optic  
Cable




Front mount & Fiber  
Optic Cable



Bullet sample without Reflector (left) and  
with Reflector (right)

## OTHER ACCESSORIES




**LCF3-SCALE**

Upgraded scale to replace older versions of LCF scales.



**LCF-COVER**

Dust Cover for LCF



**Leeds Calibration Kit**  
Magnification matching kit

## CASTERS | ANTI-VIBRATION FEET

Anti-vibration Pads, designed for under the table legs of the LCF2 and LCF3. Sold individually. Part Number: LCF3-AVPAD



Leveling Caster, Sold as set of 4



Anti-vibration feet, Sold as set of 4 for LCF or LCF3



High deflection anti-vibration feet, Sold as set of 4 for LCF3



**Leeds Precision Instruments, Inc.**  
DBA Leeds Forensic Systems  
17300 Medina Road, Suite 600  
Minneapolis, MN 55447 USA  
+1.763-546-8575 // [www.leedsmicro.com](http://www.leedsmicro.com)