



# Fluorescent Quad Lamp


## Bright, Diffuse Light

The Leeds Fluorescent Quad Lamp provides bright diffuse illumination. The **unique rotating hood** allows examiners to control the angle of light and adjust the contrast on a sample.

The Quad Lamp can be attached with an articulated arm, allowing the user to adjust and tighten joints as desired.

Or, the Quad Lamp can be attached to the microscope mounted on obedient gooseneck tubing, with stay-in-place capabilities. Adjust the tubing with one hand, and the quad lamp will remain in the desired location without needing to be locked into place.

These replacement lamps are quality control checked and sorted by their background color, to minimize the large color variance that can occur between individual lamps.




Fluorescent Quad Lamp mounted on the LCF3 Firearms and Tool Marks Comparison Microscope



Fluorescent Quad Lamp with articulated arm, attached to Focus Mount



Fluorescent Quad Lamp replacement bulbs, sold in a set of 2, on Leeds' ebay store, accessible via [www.leedsmicro.com](http://www.leedsmicro.com)



Fluorescent Quad Lamp mounted to column with flexible tube arm

Discuss your lighting needs  
with our knowledgeable staff:  
[sales@leedsmicro.com](mailto:sales@leedsmicro.com)

LEEDS PRECISION INSTRUMENTS, INC.  
17300 MEDINA ROAD, SUITE 600  
MINNEAPOLIS, MN 55447 USA  
763-546-8575  
[WWW.LEEDSMICRO.COM](http://WWW.LEEDSMICRO.COM)