





Accessories for Leeds LCF Firearms & Tool Marks Comparison Microscopes

Coaxial Illumination Kit -- an ideal solution when looking at granular breach face marks and analyzing firing pin impressions. This Kit provides glare-free illumination.



Leeds Quick Change Kit: Add this mounting kit to **easily and quickly swap out objectives** and the coaxial illuminator.



CMOS5-PKG CAMERA
Leeds CMOS5 Camera is optimized for the rounded shape of bullets and cartridge cases. It has a unique high dynamic range (HDR) mode to eliminate glare, and an extended depth of field (EDF) capability to enhance depth of focus.




Leeds Ultra-Widefield Objective


0.18x objective, large FOV at 200mm.
Continuous 1.08x to 18x parfocal zoom range.
Working Distance of 106.5mm.



SAMPLE HOLDERS



Magnetic Shell Holder




White Contrast Disk for Universal Holder



Mega Holder



Wax Platen with ball pivot



Ball Pivot Positioner with 4" diameter platen




Ball Pivot Vise Holder



Illuminator Reflector mounted on flexible tube arm.




On-Centers Holder



Reflected Light Document Holder

POST MOUNTS / SHELVING / PIN MOUNTS & EYEPIECE HOLDER

Monitor Post Mount attaches to the back of the table and can be oriented to the right, left or center of the table.




Keyboard & PC Mounts easily attach to the table leg (either right or left) and can be adjusted to desired height.




LCF3 Table Extension Kit
Extension table fits the left or right side of the LCF3 station. Includes a neoprene pad insert and our unique anti-roll-off lip.



Handy Shelf attachments, in 2 sizes



PIN MOUNT & LEEDS EYEPIECE HOLDERS



Pin Mount Holder Tray (shown here without the eyepiece holder) attaches to the Z column.




Set of Pin Mounts (Pairs)



Pin Mount Tray with a handy space for storage of your 3rd eyepiece.

LIGHTING OPTIONS

LEEDS FLUORESCENT QUAD-TUBE LIGHTS




LCF2 Fluorescent Light Cross Bar Mount, with adjustable arm ('squeaky tube')




LCF2 Fluorescent Kit mounts on Cross Bar or Saddle with Adjustable Arm. (Position light first, then tighten star knob to lock in place)



STAGE or CROSS BAR MOUNT



COLUMN MOUNT



Bifurcated Fiber Optic Light Guide with obedient gooseneck



Diffused Spot Lens



Focusing Spot Lens

QUAD-TUBE REPLACEMENT



Fluorescent Quad-tube, color-matched replacement bulbs. Sold as set of 2.


TWIN-TUBE REPLACEMENT




Twin-tube Fluorescent replacement bulbs.



Light Guide Positioner: Used to record and duplicate the distance from sample, height, angle of illumination, and rotary position. Useful for reverse lighting technique.



OTHER ACCESSORIES




LCF3-SCALE

Upgraded scale to replace older versions of LCF scales.
(Note: Not for motorized LCF)



LCF-COVER


Dust Cover for LCF



Leeds Calibration Kit:

200 division stage micrometer; reticle with 200 divisions; slide holder; flash drive with instructions & templates

X Y Z CONTROLLERS




To control the X, Y & Z stage movements on your Leeds microscope with motorized stages, choose either:

- the Human Interface Device (HIDS) (shown above) or
- the Trackball (shown below.)



CASTERS / ANTI-VIBRATION FEET



- A) Anti-vibration Pads, designed for under the table legs of LCF systems. Sold individually. P/N: LCF3-AVPAD
- B) Leveling Caster, Sold as set of 4.
- C) Anti-vibration feet, Sold as set of 4.
- D) High deflection anti-vibration feet, Sold as set of 4.