



LIGHT GUIDE POSITIONER




ACCESSORY FOR
LIGHTING REPEATABILITY



Leeds Light Guide Positioner (LGP) is an optional accessory to add to your comparison microscope configuration for **lighting repeatability**.

LGP capabilities:

- **Can record and duplicate** the distance from the sample; the height and angle of illumination; and the rotary position.
- **Reverse Lighting** -- the examiner can set up lighting from the same length, height, and angle of illumination of the sample on both sides (mirrored illumination of each other.) This feature assures that both samples are illuminated in the exact same manner.
- Ability to **rotate the illumination source** around the sample 360°.

An LGP can be used on each side of a Leeds LCF or Discovery comparison microscope.





**Light Guide Positioner for
Leeds Comparison Microscopes,
for lighting repeatability**



LEEDS PRECISION INSTRUMENTS, INC.
17300 MEDINA RD, SUITE 600
MINNEAPOLIS MN 55447 USA
+1-763-546-8575 // WWW.LEEDSMICRO.COM