

Leeds LEDFO Cube Illuminator

The **Leeds LEDFO Cube Illuminator** is an LED fiber optic light source designed for a wide variety of microscope applications.

- **Vibration-Free Operation** – Designed with no fan or moving parts.
- **Constant Color Temperature** – Maintains consistent color at all light intensities.
- **Long LED Life** – Provides an estimated **30,000 hours of operation**.
- **Cool Illumination** – Designed with no thermal component in the emitted light.
- **Versatile Compatibility** – Accommodates most manufacturers' fiber optic ring lights and light guides.
- **Low Power Consumption** – Energy-efficient design for reliable everyday operation.

