


OPTIONAL ACCESSORIES FOR FIREARMS & TOOL MARKS COMPARISON MICROSCOPES

ENHANCE **PERFORMANCE**.

IMPROVE **ACCURACY**.

SUPPORT EVERY **EXAMINATION**.



ENHANCED CLARITY

Reduce glare and
reveal more detail.



EASY INTEGRATION

Seamlessly integrates
with Leeds Comparison
Microscopes.





PRECISION OPTICS

Capture the detail that
makes the difference.

COAXIAL ILLUMINATION KIT

Ideal for specular reflection, (mirror-like) on
breach face marks and firing pin impressions.
Provides glare-free illumination for superior
detail and contrast.



LCF-CMOS5-HDR CAMERA

The Leeds CMOS5 Camera is
a high Quality Imaging Camera
featuring a fast refresh rate and
HDR capabilities to enhance image
detail when examining Bullets and
Cartridge Case.



LEEDS QUICK CHANGE KIT

Add this mounting kit to easily and quickly
swap out objectives.




ULTRA-WIDEFIELD OBJECTIVE


0.18x objective with an
Extra-wide 170 mm
Field of view, parfocal
throughout entire zoom ranges,
and 106.5 mm working distance.
(Easily Swap objectives with
Quick Change Kit)



► FIREARMS & TOOL MARKS SAMPLE HOLDERS




MAGNETIZED SHELL HOLDER



ON-CENTER HOLDER



MEGA UNIVERSAL HOLDER




BALL PIVOT WAX PLATEN



BALL PIVOT POSITIONER



**REFLECTED LIGHT
DOCUMENT HOLDER**



BALL PIVOT VISE HOLDER



BALL PIVOT POSITIONER



WAX PLATTEN HOLDER

▶ LIGHTING OPTIONS

LEEDS LED LIGHT STIX



**FLUORESCENT
QUAD-TUBE LIGHTS**



**LIGHT GUIDE
POSITIONERS**




**SIX-BUTTON KEY
PAD CONTROLLER**



ILLUMINATOR REFLECTOR




**STAGE or CROSS-BEAM
MOUNT**
articulated arm attachment
LCF3-LEDKIT-AA-R2



**STAGE or CROSS-BEAM
MOUNT**
flexible tube arm,
cross-beam mount
LCF3-LEDKIT-CB



**QUAD-TUBE
REPLACEMENT**



**TWIN-TUBE
REPLACEMENT**



FIBER OPTIC LIGHTING

▶ POST MOUNTS, SHELVING & TABLE EXTENSIONS



**MONITOR POST
MOUNT**



**KEYBOARD &
PC MOUNTS**



**LCF 16" SHELF
FOR LCF3**




**LCF3 TABLE
EXTENSION KIT**

▶ CASTERS & ANTI-VIBRATION FEET



LCF3 Anti-Vibration Pads

Sold in set of 4, reduce
vibration under LCF3 legs



LCF3-LVCSTR
Leveling Caster, Set of 4



LCF3-VIBR
Anti-vibration feet, set
of 4 for LCF or LCF3



LCF3-VIBR-.5
High deflection anti-vibration
feet, set of 4 for LCF3