Fluorescent Quad Lamp
BRIGHT, DIFFUSE LIGHT

40th Anniversary
LEEDS
FORENSIC SYSTEMS
Your Forensic Imaging Source

Fluorescent Quad Lamp

The Leeds Fluorescent Quad Lamp provides bright, diffuse illumination for examining a wide variety of samples. Its unique rotating hood allows examiners to control the angle of light and adjust sample contrast, helping reveal surface details during examination.


The Quad Lamp offers two flexible mounting options to accommodate different positioning needs. An articulated arm allows users to adjust the lamp's position and tighten the joints to hold it at the desired angle.

Alternatively, the lamp can be mounted to the microscope using stay-in-place goose-neck tubing. This option allows easy, one-handed adjustment and keeps the lamp in the desired location without requiring joints to be tightened or locked.





Fluorescent Quad Lamp replacement bulbs, sold in a set of 2, on Leeds' ebay store, accessible via **www.leedsmicro.com**.







Fluorescent Quad Lamp mounted on the LCF3 Firearms and Tool Marks Comparison Microscope.



Fluorescent Quad Lamp with articulated arm, attached to Focus Mount.

www.leedsmicro.com **Phone: 763-546-8575** **Email: sales@leedsmicro.com**